

FS-B Series

FS-B series are excellent products which are made by high technology of 0.01g accuracy, has feature of standard wind shield which is not use in same level product, response quickly, high sensitivity and stability.



CHARACTERISTIC

- > LCD back Light
- > Response time in 2 seconds
- > Support AC-DC
- > Selectable weighing units and count mode
- > 10 unit available: g, ct, lb, oz, ozt, GN, dwt, N,t, T/A/R
- > Fit with RS232 Line
- > Fit with windshield
- * Fit with inner rechargeable DC battery(Optional)
- * Fit with extra display(Optional)
- * Built in Hook(Optional)

SPECIFICATION

Type	Metage Range(g)	Reading Precision(g)	Repeatability (g)	Linearity(g)	Ambient temperature(°C)	Scale Pan Size(mm)	Shape Size(LxWxH)(mm)
FS-B100	110	0.01	±0.01	±0.02	5-35	Φ128	270x190x265
FS-B200	210	0.01	±0.01	±0.02	5-35	Φ128	270x190x265
FS-B300	310	0.01	±0.01	±0.02	5-35	Φ128	270x190x265
FS-B500	510	0.01	±0.01	±0.02	5-35	Φ128	270x190x265
FS-B600	610	0.01	±0.01	±0.02	5-35	Φ128	270x190x265
FS-B1000	1100	0.01	±0.01	±0.02	5-35	Φ138	270x190x265
FS-B2000	2100	0.01	±0.02	±0.02	5-35	140x160mm	270x190x265